

SUREFIRE® INFUSION SYSTEM (SIS) SPECIFICATIONS

Product Types	Ordering Code	Inner Diameter	Vessel Size	Usable Length	Proximal Outer Diameter	Distal Tip Outer Diameter	Dead Space	Max Contrast Flow Rate (@1200psi)	Bead Size Compatibility
---------------	---------------	----------------	-------------	---------------	-------------------------	---------------------------	------------	-----------------------------------	-------------------------


120 cm CATHETERS

021	SPS-21120-35	0.021 in	1.5 - 3.5 mm	120 cm	3.2 F	2.8 F	0.37 mL	2.6 mL/sec	Hydrogel ≤ 500 µm Glass ≤ 110 µm
025M	SPS-25120-40	0.025 in	2.0 - 4.0 mm	120 cm	3.7 F	3.4 F	0.52 mL	3.6 mL/sec	Hydrogel ≤ 500 µm Glass ≤ 150 µm
025L	SPS-25120-60	0.025 in	4.0 - 6.0 mm	120 cm	3.7 F	3.4 F	0.52 mL	3.6 mL/sec	Hydrogel ≤ 500 µm Glass ≤ 150 µm

150 cm CATHETERS

021	SPS-21150-35	0.021 in	1.5 - 3.5 mm	150 cm	3.2 F	2.8 F	0.44 mL	2.1 mL/sec	Hydrogel ≤ 500 µm Glass ≤ 110 µm
025M	SPS-25150-40	0.025 in	2.0 - 4.0 mm	150 cm	3.7 F	3.4 F	0.61 mL	3.3 mL/sec	Hydrogel ≤ 500 µm Glass ≤ 150 µm
025L	SPS-25150-60	0.025 in	4.0 - 6.0 mm	150 cm	3.7 F	3.4 F	0.61 mL	3.3 mL/sec	Hydrogel ≤ 500 µm Glass ≤ 150 µm

All catheters listed are incompatible with glue or alcohol.



ACCESSORY: BASE CATHETERS

Speak with a representative for a complete list of compatible shapes or other base catheters.

Manufacturer	Brand	Shape	Size	Inner Diameter	Usable Length
Surefire Medical	Flexion	Axis, SIM1	5F	0.054 in	80 cm
Surefire Medical	Surefire®	Axis, Cobra1, SIM1	5F	0.056 in	65 cm

Compatible with SIS 021, 025M & 025L 150 cm Catheters

DePuy Synthes/ Codman	Envoy®	Multipurpose C and D, Modified Cerebral and Headhunter 1	5F	0.056 in	90 cm, 100 cm
Medtronic	Launcher®	JR 4.0, MB1a and PAPA	5F	0.058 in	90 cm, 100 cm, 110 cm
Merit Medical	ConcierGE®	All Ultimate™ Curves and Multipurpose, Hockey Stick and Judkins Right 3.5	5F	0.057 in	100 cm
Terumo	Heartrail® III	TIG4.0, Ikari Right and Ikari Left	5F	0.058 in	100 cm

CONTACT US:

trisaluslifesci.com
info@ trisaluslifesci.com

Toll Free +1.888.321.5212
Office +1.303.426.1222

Angiographic catheters with inner diameters from 0.035-0.038 in are not compatible.